



200/75 PRIME M

Cod. PR-OF900-200-MT3

Oil free reciprocating air compressors with a compressed air FAD of 23.3 CFM @ 72.5 PSI pressure. Receiver capacity 200 L, with internal surface treatment to protect against corrosion.

Complete with membrane dryer to ensure maximum air purity. This compressor is available also in [145 PSI VERSION](#), [visit the appropriate section!](#)













Maximum simultaneous use: 9-12




Dental Chairs: 12

Medical Device CE 0425

Cod.	Air Dryer	 Air capacity	 FAD Fad @ 72.5 PSI	 Motor	 Voltage Hz	 Receiver	 Noise Level	 Duty Cycle	 max P.	 Size	 net Weight
200/75 TANDEM PRIME M PR-OF900-200-T3	Yes	cfm 38.12	cfm 23.3	HP 8.6	400V 50Hz	GAL 52.84	db(A) 72	70 %	PSI 116	inch 57.1x32.3x35.4	Lb 403.5
200/75 TANDEM PRIME M UL PR-OF900-200-T3-UL	Yes	cfm 38.12	cfm 23.3	HP 8.6	400V 50Hz	GAL 52.84	db(A) 72	70 %	PSI 116	inch 57.1x32.3x35.4	Lb 403.5

Air delivery performance



Maintenance kits

Kit.	Description	Type	Hours
TA-FL-0263-K-SCP	<p>TANDEM PRIME 3 + 3 + 3 M kit</p> <p>Ordinary maintenance kit for the PURE AIR TANDEM PRIME 3+3+3 M range (200/75 PRIME M) with air dryer</p> <ul style="list-style-type: none"> • 9 suction filters, code: TA-FL-058-C • 2 non-return valve pads, cod: CO-PV-001 • 3 5 micron filter cartridges, code: TA-FL-001-C-SM • 3 0.01 micron filter cartridges, code: TA-FL-024-C2-SM • 6 o-rings for filters, cod. TA-FL-001-OR • 3 radiator filter cartridges, code: TA-FL-002-C • 3 0.01 micron filter cartridges, code: TA-FL-076-C • 3 AMG150 water separator cartridges, code: TA-FL-074-C • 6 o-rings for filters, cod. TA-FL-075-OR <p>Replace the filtration cartridges every 2 years.</p>	Ordinary	1.500h/2y
PR-GC-035-K-SCP	<p>TANDEM PRIME 3+3+3 Kit</p> <p>The extraordinary maintenance spare parts kit for the PRIME 3+3+3 versions (200/75 PRIME).</p> <p>Comprehends:</p> <ul style="list-style-type: none"> • 9 suction filters, code: TA-FL-058-C • 18 O-rings, code: MA-PC018 • 9 piston rings , code: MA-MG00001 • 9 piston guide rings, code: MA-MG00003 • 9 expansion rings for the segments, code: MA-MG00024 • 9 complete valve plates, code: MA-MG-02002P-C-P • 3 non-return valve pads, cod: CO-PV-001 <p>The recommended maintenance intervals for this compressor are 4000 hours/4 years.</p>	Extraordinary	4.000h/4y